

Comprehensive Product Range

Totalling hundreds of models of Sharps Containers.

New Designs

More than 10 types of new designs every year.

Quick Shipment

Minimum delivery time is 3 days.

High Performance & Low Cost

Let you be very satisfied with the products.

Service Guarantee

Response within 24 hours for all questions.



Available Colour:

Contents

Company Profile	3
Product Application	4
K Series Sharps Container	
K0.2	5
F Series Sharps Container	
Introduction	6
F02 F03	7
F04 F05A	8
F05B F06	9
F07 F10	10
F15 F23	11
R Series Sharps Container	
Introduction	12
R0.7 R1	13
R1.5 R2.4	14
R3A R3B	15
R6 R8	16
Y Series Sharps Container	
Introduction	17
Y01 Y03	18
Y04 Y06A	19
Y06B	20
T Series Sharps Container	
Introduction	21
T1A T1B	22
T2 T3	23
T4.6A T4.6B	24
T4.6C T5	25
T8	26
S Series Sharps Container	
Introduction	27
S1.7 S2	28
S7 SPS	29
Big size Sharps Container	30
Respiratory exerciser	30





Dailymag is China's leading manufacturer and exporter of sharps containers.

Utilizing advanced production equipment and a rigid quality control system, we provide five different series of sharps containers, totaling hundreds of models. We take great pride in the fact that our products protect both workers and the general public from unnecessary injury or possible infectious disease transmission.

We are keenly focused on applying our advanced engineering and production technology to provide both standardized and custom design solutions that lower costs and supply comprehensive sharps container solutions to our worldwide customers. We promise and deliver to each client a comprehensive product range, new designs, quick shipment, high performance, and low cost.

Dailymag is ISO 9001:2008 Certified, and our products are BS7320, UN and FDA approved. We have more than 100 patents, including more than 20 for sharps containers.

For Dailymag, sustainability is a word that can be used to describe our business relationships, as we value our clients and consistently deliver value to them. To that end we are excited about the prospect of establishing and fostering a mutually beneficial long-term relationship with you and your company, and look forward to that beginning.



Dailymag Market



Application of Medical Sharps Container

Medical sharps containers are for properly dispose of medical waste which includes but not limited to, discarded medical needles, scalpels, IV catheters, non-liquid human tissue, body parts, human blood, human blood products, laboratory waste which contain human disease-causing agents, and body fluids.

We use sharps container to ensure all the waste are properly segregated from people to prevent infection and harm in hospitals, clinics, labs and also home care.



K Series Pocket Sharps Container

- Smart size
- Leak proof and puncture resistant
- Appropriate sliding opening
- Temporary closure & definite closure
- Needle removal cap
- Befitting for the pocket
- Very convenient for outdoor or travel use



Model: DMS-K0.2 (Capacity: 0.2L)

Product size (cm)	pcs/case	Carton size (cm)
7.6x3.8x10	200	30x35x42



F Series

- Most popular quadrate shape
- Leak proof and puncture resistant
- Wide and clear opening for easy disposal
- Counterbalanced door for hands-free disposal
- Ergonomically-designed handle
- Temporary closure & definite closure
- Practical rotating/sliding temporary lid
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment



www.sharp-container.com

F Series

Model: DMS-F02 (Capacity: 2.0L)

Product size (cm)	pcs/cbm	Carton size (cm)
20x16x12	30	41x34x23



F Series

Model: DMS-F04 (Capacity: 4.0L)

Product size (cm)	pcs/cbm	Carton size (cm)
28.5x20x14.5	24	51x42x25



Model: DMS-F03 (Capacity: 3.0L)

Product size (cm)	pcs/cbm	Carton size (cm)
20x16x13.5	30	41x34x31.5



Model: DMS-F05A (Capacity: 5.0L)

Product size (cm)	pcs/cbm	Carton size (cm)
28.5x20x16.5	24	51x42x29.5



F Series

Model: DMS-F05B (Capacity: 5.0L)

Product size (cm)	pc/cbm	Carton size (cm)
24.5x20x16	24	51x42x27



Model: DMS-F06 (Capacity: 6.2L)

Product size (cm)	pc/cbm	Carton size (cm)
24.5x20x18.5	24	51x42x31



F Series

Model: DMS-F07 (Capacity: 7.0L)

Product size (cm)	pc/cbm	Carton size (cm)
25x20x23	24	51x42x37



Model: DMS-F10 (Capacity: 10.0L)

Product size (cm)	pc/cbm	Carton size (cm)
30x20x26	12	62x34x38



F Series



Model: DMS-F15 (Capacity: 15.0L)

Product size (cm)	PCS/CM	Carton size (cm)
36x22x36	20	49x37x65



Model: DMS-F23 (Capacity: 23.0L)

Product size (cm)	PCS/CM	Carton size (cm)
37x24x39	20	49x37x62



R Series

- New-developed round shape
- Leak proof and puncture resistant
- Funnel-designed and built-in needle removal cap for easy and safe disposal
- Various indent design for easy needle removal
- Temporary closure & definite closure
- Convenient extra lid design for temporary and final closure
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment

Available Colour: 

R Series

Model: DMS-R0.7 (Capacity: 0.7L)

Product size (cm)	pcs/ctn	Carton size (cm)
9x9x7.5x13.5	100	39x33x99



Model: DMS-R1 (Capacity: 1.0L)

Product size (cm)	pcs/ctn	Carton size (cm)
11x8.5x15	100	39x33x58



R Series

Model: DMS-R1.5 (Capacity: 1.5L)

Product size (cm)	pcs/ctn	Carton size (cm)
14x11x15	80	37x32x54



Model: DMS-R2.4 (Capacity: 2.4L)

Product size (cm)	pcs/ctn	Carton size (cm)
14x11.5x20	40	48x32x50



R Series



Model: DMS-R3A (Capacity: 3.0L)

Product size (cm)	pa3/cm	Carton size (cm)
14x12x24	30	40x32x65



R Series



Model: DMS-R6 (Capacity: 6.0L)

Product size (cm)	pa3/cm	Carton size (cm)
22.5x20x19	24	44x40x62.5



Model: DMS-R3B (Capacity: 3.0L)

Product size (cm)	pa3/cm	Carton size (cm)
14x12x24	30	40x32x65



Model: DMS-R8 (Capacity: 8.0L)

Product size (cm)	pa3/cm	Carton size (cm)
22.5x18x26	24	53x42x65



Y Series

- Traditional-designed round shape
- Leak proof and puncture resistant
- Wide opening for easy disposal
- Temporary closure & definite closure
- Large and practical rotating/sliding temporary lid
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment

Available Colour: 

Y Series

Model: DMS-Y01 (Capacity: 0.8L)

Product size (cm)	Pcs/ctn	Carton size (cm)
126x115	90	1843x915



Model: DMS-Y03 (Capacity: 2.8L)

Product size (cm)	Pcs/ctn	Carton size (cm)
186x147.5	30	1863x915



Y Series




Model: DMS-Y04 (Capacity: 3.6L)

Product size (cm)	pcs/cdn	Carton size (cm)
18.3x19x19	30	32x32x32



Model: DMS-Y06A (Capacity: 6.0L)

Product size (cm)	pcs/cdn	Carton size (cm)
22x19x20	30	44x44x44.3



Y Series



Model: DMS-Y06B (Capacity: 6.0L)

Product size (cm)	pcs/cdn	Carton size (cm)
22x19x20	30	44x44x44.3



T Series

- Novel quadrate shape
- Leak proof and puncture resistant
- Wide, clear and practical opening for easy and safe disposal
- Counterbalanced door for hands-free disposal
- Built-in needle removal
- Ergonomically-designed handle on the side/lid
- Temporary closure & definite closure
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment



Available colour: 




T Series

Model: DMS-T1A (Capacity: 1.0L)

Product size (cm)	pcs/ctn	Carton size (cm)
10.5x8x18	60	37x32x41



Model: DMS-T1B (Capacity: 1.0L)

Product size (cm)	pcs/ctn	Carton size (cm)
11x10x15.5	60	37x32x41



T Series



Model: DMS-T2 (Capacity: 2.0L)

Product size (cm)	pl/ctn	Carton size (cm)
79x11x16.5	30	40x32x74x3



Model: DMS-T3 (Capacity: 3.0L)

Product size (cm)	pl/ctn	Carton size (cm)
16x13x18	54	47x35x40

T Series



Model: DMS-T4.6A (Capacity: 4.6L)

Product size (cm)	pl/ctn	Carton size (cm)
27x11.5x28	20	58x37x29



Model: DMS-T4.6B (Capacity: 4.6L)

Product size (cm)	pl/ctn	Carton size (cm)
27x11.5x27	20	58x35x38



S Series

- Special and square shape
- Newly designed
- Leak proof and puncture resistant
- Solid and stable stand
- 2L and 5L models are perfect for long size needle disposal
- Cardboard sharps box for special needs
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment

Available Colours: 

S Series

Model: DMS-S1.7 (Capacity: 1.7L)

Product size (cm)	pcs/ctn	Carton size (cm)
10x10x21.5	40	54x44x46.5



Model: DMS-S2 (Capacity: 2.0L)

Product size (cm)	pcs/ctn	Carton size (cm)
13x12x25.5	40	40x33x49



S Series



Model: DMS-S7 (Capacity: 7.0L)

Product size (cm)	pc/s/ctn	Carton size (cm)
20x20x26.5	24	44x42x45



Model: DMS-SP5 (Capacity: 5.0L)

Product size (cm)	pc/s/ctn	Carton size (cm)
15x12.5x29	60	47x42x28.5



Big Size Sharps Container

• Available capacity: 25L-240L



Respiratory exerciser

- Help patients perform and monitor their own post-surgical breathing exercisers
- One-way valve ensures patients inhale rather than exhale into the unit
- Highly visible pistons and universal graphics indicate correct inspiratory flow rate

Model: DMS-A1
(Capacity: 600cc/900cc/1200cc per sec)

Product size (cm)	pc/ctn	Carton size (cm)
16x7x14	50	30x26x29.5


www.sharp-container.com

Dailymag[®]
Your Best Source in Sharps Containers