






AGS-MEDICAL




Ph: 505-550-6501 & 505-565-5102, Fx: 505-814-5778





Em: sales@agsmedical.com, Web: http://www.agsmedical.com






Item No.	Picture	FOB	Specification	set	Details
JQ-FM-5 Electric weighing bed(Five function)		\$	L2170XW1080XH500/850mm In bed plate:L1950XW910mm 	Overall lifting 500-850mm,back plate 65°±5°,leg plate 45°±5°,Head and Foot tilt 13°±3	1.The whole bed can weight,with error controlled in 1-3% 4 pcs of european style guardrail,upside-down structure, 3Imported Germany "TENTE"central control castors with lock,safety and mute.4.Bed board adopts HPL board,X-ray photograph for whole body is available 5.With anti-collision wheel and IV pole(optional) the drainage hook at two sides
JQ-FM-4 Electric bed(Five function)		\$	L2170XW1080XH440/800mm In bed plate:L1950XW910mm 	Overall lifting 440-800mm,back plate 70°±5°,leg plate 45°±5°,tilt:16°±3°	4 pcs of european style guardrail,upside-down structure, luxury central control castors,with lock,stright going,free moving braking device, bed board adopts HPL board,X-ray photograph for whole body is available, 2 step adjustment at shank section
JQ-FM-3 Electric bed (Four function)		\$	L2170XW1080XH400/750mm In bed plate:L1950XW910mm 	Overall lifting 400-750mm,back plate 70°±5°,leg plate 45°±5°	4 pcs of european style guardrail,upside-down structure, luxury central control castors,with lock,stright going,free moving braking device,
JQ-FN-2 Manual bed (three function)		\$	 L2200XW1020XH420/760mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:420-760mm	High quality ABS engineering plastic head/foot boards,easily detachable and stable 2.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust. 3.ABS folded cranks Luxurious central locak castors with diameter 125mm

Item No.	Picture	FOB	Specification	set	Details
JQ-FN-1 Manual bed (two function)			 Size:L2200XW1020XH500mm	back rest:0°-75° knee rest:0°-45°	High quality ABS engineering plastic head/foot boards,easily detachable and stable 2.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust. 3.ABS folded cranks Luxurious central locak castors with diameter 125mm
JQ-FE Series (Electric & Manual)					
JQ-FE-1A electric bed (four function)			 L2175XW1090XH420/720mm	Overall lifting 420-720mm, back plate 0°-75°, leg plate 0°-45°	1.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust.2.High quality ABS engineering plastic head/foot boards,easily detachable and stable 3.Safety voltage motor,noiseless
JQ-FE-1 electric bed (three function)			 L2175XW1090XH420/720mm	Overall lifting 420-720mm, back plate 0°-75°, leg plate 0°-45°	4.Equipped with diameter 125mm luxury mute central controlled castors 5. IV pole (optional) the drainage hook at two sides of the bed
JQ-FE-2 Manual bed (three function)			 L2175XW1090XH420/720mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:420-720mm	1.ABS wall mounted bed head and end,easy for mounting and dismounting(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails3.Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller


Item No.	Picture	FOB	Specification	set	Details
JQ-FE-3 Manual bed(two function)			 Size:L2175XW1090XH500mm	Manual control back rest:0°-75°, knee rest:0°-45°	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails3.Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand
JQ-FA Series (Electric & Manual)					
JQ-FA-1 electric bed(three function)			 L2200XW980XH380/600mm	Overall lifting 380-600mm, back plate 70°±5°, leg plate 40°±5°	1.the handle is made by ABS injection forming of high strength at a time 2.concave design with smooth punched holes,easy for cleaning and sorting 3.ABS lifting protectng rails,in the same color with bed head and end,luxurious and beautiful 4.Luxurious central locak castors with diameter 125mm
JQ-FA-2 Manual bed(three function)			 L2200XW980XH380/600mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:380-600mm	1.the handle is made by ABS injection forming of high strength at a time 2.concave design with smooth punched holes,easy for cleaning and sorting 3.ABS lifting protectng rails,in the same color with bed head and end,luxurious and beautiful 4.Luxurious central locak castors with diameter 125mm
JQ-FA-3 Manual bed(two function)			 L2200XW980XH500mm	Manual control back rest:0°-75°, knee rest:0°-45	1.the handle is made by ABS injection forming of high strength at a time 2.concave design with smooth punched holes,easy for cleaning and sorting 3.ABS lifting protectng rails,in the same color with bed head and end,luxurious and beautiful 4.Luxurious central locak castors with diameter 125mm
JQ-FA-4 Manual bed(two function)			 L2200XW980XH500mm	Manual control back rest:0°-75°, knee rest:0°-45	1.the handle is made by ABS injection forming of high strength at a time 2.concave design with smooth punched holes,easy for cleaning and sorting 3.ABS lifting protectng rails,in the same color with bed head and end,luxurious and beautiful 4.Four foot brake castor




Item No.	Picture	FOB	Specification	set	Details
JQ-FA-5 electric bed(three function)			 Size:L2200XW960XH420/720m	Overall lifting 420-720mm, back plate 70°±5°, leg plate 40°±5°	1.Aluminum alloy guardrails,with finger protecting design 2. Luxury central control castors,brake the bed with one step. 3.Bed panel is made of one-step molding cold rolled steel,excellent antu-slipping performance 4.Load capacity over 240kg
JQ-FA-6 Manual bed(three function)			 L2200XW960XH440/790mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:440-790mm	1.High quality ABS engineering plastic head/foot boards,easily detachable and stable. 2.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust. 3.ABS folded cranks 4.Luxurious central locak castors with diameter 125mm
JQ-FA-7 Manual bed(two function)			 L2200XW960XH520mm	Manual control back rest:0°-75°, knee rest:0°-45	1.High quality ABS engineering plastic head/foot boards,easily detachable and stable. 2.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust. 3.ABS folded cranks 4.Luxurious central locak castors with diameter 125mm
JQ-FA-8 Manual bed(two function)			 L2200XW960XH520mm	Manual control back rest:0°-75°, knee rest:0°-45	1.High quality ABS engineering plastic head/foot boards,easily detachable and stable. 2.Bed frame welded by profile steel,bed platform is treated with advanced eletrostatic plastic-spray which is antiaging and anti-rust. 3.ABS folded cranks 4.Four foot brake castor
JQ-FB Series (Electric & Manual)					


Item No.	Picture	FOB	Specification	set	Details
JQ-FB-1 Electric bed(three function)			 L2200XW950XH380/600mm	Overall lifting 380-600mm, back plate 70°±5°, leg plate 40°±5°	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails 3.Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FB-2 Manual bed(three function)			 L2200XW950XH380/600mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:380-600mm	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails 3.Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FB-3 Manual bed(two function)			 L2200XW950XH530mm	Manual control back rest:0°-75°, knee rest:0°-45	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails 3.Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FB-4 Manual bed(two function)			 L2200XW950XH530mm	Manual control back rest:0°-75°, knee rest:0°-45	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy protecting rails,available for concealing into the bed for optionally selecting table board adaptive to protecting rails 3. Four foot brake castor 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FC Series (Electric & Manual)					




Item No.	Picture	FOB	Specification	set	Details
JQ-FC-1 electric bed(three function)			 L2150XW950XH380/600mm	Overall lifting 380-600mm, back plate 70°±5°, leg plate 40°±5°	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy anti-pinch protecting rails,available for concealing into the bed,table board adaptive to protecting rails (optional) 3. Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FC-2 Manual bed(three function)			 L2150XW950XH380/600mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:380-600mm	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy anti-pinch protecting rails,available for concealing into the bed,table board adaptive to protecting rails (optional) 3. Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FC-3 Manual bed(two function)			 L2200XW950XH530mm	Manual control back rest:0°-75°, knee rest:0°-45	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy anti-pinch protecting rails,available for concealing into the bed,table board adaptive to protecting rails (optional) 3. Central-control brake wheel,with higher stability 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller
JQ-FC-4 Manual bed(two function)			 L2150XW950XH530mm	Manual control back rest:0°-75°, knee rest:0°-45	1.ABS wall mounted bed head and end,easy for mounting and dismantling(with safety pin) 2.Aluminum alloy anti-pinch protecting rails,available for concealing into the bed,table board adaptive to protecting rails (optional) 3. Four foot brake castor 4.Be equipped with Denmark safety-voltage motor(LINAK) omnibearing and adjustable hand controller



JQ-FP Series (Electric & Manual)


Item No.	Picture	FOB	Specification	set	Details
JQ-FP-7 Electric bed(seven function)			 L2240XW1100XH630/1030mm	Overall lifting 630-1030mm,back plate 65°,leg plate 35°,Trendelengurg angle $\geq 12^\circ$ Reverse trendelengurg angle $\geq 15^\circ$,lateral tilt $\geq 12^\circ$	1 pair of european style bed head; ABS damping lifting guardrail; Soft connection bed penal; 5 pcs of Motor; 4 control panels inside and outside guardrails; 4 Imported German "TENTE" medical ultra silent central castor; Multiple position electric adjusting function,enabling heart displacement,shock,CPR and examination positions;overall waterproof level is
JQ-FP-5 Electric bed(Five function)			 L2240XW1100XH630/1030mm	Overall lifting 630-1030mm,back plate 65°,leg plate 35°,Trendelengurg angle $\geq 12^\circ$ Reverse trendelengurg angle $\geq 15^\circ$	1 pair of european style bed head; ABS damping lifting guardrail; Soft connection bed penal; 4 pcs of Motor; 4 Imported German "TENTE" medical ultra silent central castor;
JQ-FP-5A Electric bed(Five function)			 L2200XW1100XH425/760mm	Overall lifting 425-760mm,back plate $70^\circ \pm 5^\circ$,leg plate $45^\circ \pm 5^\circ$,tilt: $16^\circ \pm 3^\circ$	
JQ-FP-4 Electric bed(four function)			 L2200XW1020XH420/780mm	Overall lifting 420-780mm,back plate $70^\circ \pm 5^\circ$,leg plate $45^\circ \pm 5^\circ$,head amd foot tilt: $\geq 12^\circ$	
JQ-FP-3 Electric bed(three function)			 L2200XW1100XH425/760mm	Overall lifting 425-760mm, back plate $70^\circ \pm 5^\circ$, leg plate $40^\circ \pm 5^\circ$	4 pcs of european style guardrail,upside-down structure;Luxury central control castors,brake the bed with one step;bed panel is made of one-step molding cold rolled steel,with excellent performance in anti-slipping and air permeability

Item No.	Picture	FOB	Specification	set	Details
JQ-FP-2 Manual bed(three function)			 L2200XW1100XH425/760mm	Overall lifting 425-760mm, back plate 70°±5°, leg plate 40°±5°	
JQ-FP-1 Manual bed(Two function)			 L2200XW1100XH500mm	Manual control back rest:0°-75°, knee rest:0°-45	
JQ-FP-6 Manual bed(single function)			 L2200XW1100XH500mm	Manual control back rest:0°-75°,	ABS headboards,Aluminum alloy protecting rails,.Central-control brake wheel
JQ-FP-7 Manual bed(single function)			 L2200XW1100XH500mm	Manual control back rest:0°-75°,	ABS headboards,Aluminum alloy protecting rails,Four foot brake castor
JQ-FP-8 Manual bed(Two function)			 L2200XW1100XH500mm	Manual control back rest:0°-75°, knee rest:0°-45	ABS headboards,Aluminum alloy protecting rails,.Central-control brake wheel
JQ-FP-9 Manual bed(Two function)			 L2200XW1100XH500mm	Manual control back rest:0°-75°, knee rest:0°-45	ABS headboards,Aluminum alloy protecting rails,Four foot brake castor




Item No.	Picture	FOB	Specification	set	Details
JQ-FO Series (Electric & Manual)					
JQ-FO-3 Electric bed(three function)			 2280*1120*420/780mm	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:420-780mm	
JQ-FO-2 Manual bed(three function)			 2280*1120*420/780	Manual control back rest:0°-75°, knee rest:0°-45°,high/low range:420-780mm	4 pcs of european style guardrail,upside-down structure;Luxury central control castors,brake the bed with one step;bed panel is made of one-step molding cold rolled steel,with excellent performance in anti-slipping and air permeability
JQ-FO-1 Manual bed(Two function)			 2280*1120*500	Manual control back rest:0°-75°, knee rest:0°-45	
JQ-FO-6 Manual bed(single function)			 2280*1120*500	Manual control back rest:0°-75°,	ABS headboards,Aluminum alloy protecting rails, Central-control brake wheel
JQ-FO-7 Manual bed(single function)			 2280*1120*500	Manual control back rest:0°-75°,	ABS headboards,Aluminum alloy protecting rails, Four foot brake castor

Item No.	Picture	FOB	Specification	set	Details
JQ-FO-8 Manual bed(Two function)			 2280*1120*500	Manual control back rest:0°-75°, knee rest:0°-45	ABS headboards,Aluminum alloy protecting rails,. Central-control brake wheel
JQ-FO-9 Manual bed(Two function)			 2280*1120*500	Manual control back rest:0°-75°, knee rest:0°-45	ABS headboards,Aluminum alloy protecting rails, Four foot brake castor
Multi-function series beds					
JQ-FG-5 Electric Seven-function Nursing Bed			 L2200XW1110XH630/1030mm	Overall lifting 630-1030mm,back plate 40°±2°,leg plate 13°±2°,Trendelengurg angle 17°±2° Reverse trendelengurg angle 25°±2°,lateral tilt 12°±2°	1 pair of european style bed head; ABS damping lifting guardrail; Soft connection bed penal; 5 pcs of Motor; 4 control panels inside and outside guardrails; 4 Imported German "TENTE" medical ultra silent central castor; Multiple position electric adjusting function, enabling heart displacement,shock,CPR and examination positions;overall waterproof level is IPX4;with weighing function;
JQ-MJ family bed(manual)			 2280X1120X500MM	manual contral back plate 65°, leg plate 35°	wooden headboards,easy to remove; steel mesh bedboard;Aluminum alloy protecting rails

Item No.	Picture	FOB	Specification	set	Details
JQ-FJ family bed(electric)			 2114X1050X650MM	electric control back plate 65°, leg plate 35°	wooden headboards,easy to removable;steel mesh bedboard;Danish safety voltage full range adjustable hand controls, Aluminum alloy protecting rails, four wheels with brake
JQ-AJ family bed(electric)			 2150X960X250/600mm	Overall lifting 250- 600mm, back plate 65°, leg plate 35°	wooden headboards,easy to removable;steel mesh bedboard; steel insert guardrail;Danish safety voltage motor(LINAK);full range adjustable hand controls
JQ-FH-12 Electric 12-function nursing bed			 Size:L2200XW1010XH400mm/750mm In bed plate: L1950XW830mm	Overall lifting 400- 750mm,back plate 65°±5°,leg plate 45°±5°,Head and Foot tilt 13°±3	Linear manipulator central control castor suction and suction hose Intelligent cleaning care host Anti-collision wheel include mattress and toilet system
JQ-FS-1 Electric Turn over nursing bed			L2250XW1340XH500mm		High quality ABS engineering plastic head/foot boards beautiful and generous,which is easily detachable and stable 24v DC motor, ultra-low noise, waterproof and dustproof Hidden type central control flexible turn over system,apply to ICU and Burned patients
Folding bed series					
JQ-Z-01			2000*900*450mm	foldable from middle line	Jet molding steel folding bed

Item No.	Picture	FOB	Specification	set	Details
JQ-Z-02			2000*900*450mm	foldable from middle line	stainless steel folding bed
JQ-Z-03			1900*770*470mm (1240X770X350mm)	foldable from middle line	Jet molding steel folding bed,include 4 castors with brake
JQ-FZ(folding bed)			2100*900*500mm	foldable from middle line	ABS headboard,include 4 castors with brake
JQ-A series manual type					
JQ-A24			2150*900*500mm	Flat bed	ABS hanging head,cold-rolled steel plate plat bed
JQ-A27			2150*900*500mm	Flat bed	ABS hanging head strip style flat bed

Item No.	Picture	FOB	Specification	set	Details
JQ-A33			2030*900*500mm	Flat bed	stainless steel head strip plat bed
JQ-A36			2030*900*500mm	Flat bed	carbon steel with spraying plat bed
JQ-A23			2150*900*500mm	Manual control back rest:0°-75°, side rails and castors are optional	ABS hanging head,cold-rolled steel plate single shake bed
JQ-A23		include four castors with brake,include side rails			
JQ-A26			2150*900*500mm	Manual control back rest:0°-75°, side rails and castors are optional	ABS hanging head strip style single shake bed
JQ-A32			2030*900*500mm	Manual control back rest:0°-75°	stainless steel head strip type single shake bed
JQ-A35			2030*900*500mm	Manual control back rest:0°-75°	Carbon steel with spraying single shake bed

Item No.	Picture	FOB	Specification	set	Details
JQ-A25			 2150*900*500mm	Manual control back rest:0°-75°, knee rest:0°-45°, side rails and castors are optional	ABS hanging head strip style double shake bed
JQ-A31			 2030*900*500mm	Manual control back rest:0°-75°, knee rest:0°-45°	stainless steel head strip type double shake bed
JQ-A34			 2030*900*500mm	Manual control back rest:0°-75°, knee rest:0°-45°	Carbon steel with spraying double shake bed
JQ-A37-1			 L2020XW900XH500/550mm	Manual control back rest:0°-75°, knee rest:0°-45°, high/low range:500-550mm	ABS head multi-function orthopedic fraction bed,The head and foot board made of ABS
JQ-A37-3			 L2020XW900XH500/550mm	Manual control back rest:0°-75°, tilt:5°±2° knee rest:0°-45°, high/low range:500-550mm	Stainless steel multi-function orthopedic fraction bed,The whole bed made of 304 stainless steel

Remarks: the bed packed in wooden case,A series bed packed in carton

Item No.	Picture	FOB	Specification	set	Details
----------	---------	-----	---------------	-----	---------

Warranty 1 year